

AUTOMATIC AUTOCLAVE

- Fully Automatic Vertical Steam Sterilizer
- Automated Cycle Sterilization
- Automatic Water Supply Water Level Sensor & Sub-Bucket
- Made With Permanent Steel For Internal Chamber, Door & Door Jig For Long Lasting Use



PT PURNAMA LABORATORY

SPECIFICATION

Type	Fully Automatic type			
Chamber Type	Round Vertical type			
Capacity	50 Liter	60 Liter	80 Liter	100 Liter
Temp. Range	Ambient to 130°C			
Temp. Controller	LED Microprocessor programmable control			
Time range	99 Hours 59 min Digital Display			
Temp. Sensor	Pt-100Ω (Exclusively for sterilizing, Auto tuning)			
Accuracy	±0.5°C			
Material(in) (out)	Stainless steel SS41 steel plate and powder coating			
Door system(out)	one touch handle type (single screw door)			
Basket	Stainless steel (#304) manufactured for sterilizing 2ea			
Pressure range	0~3kg/cm ² (Analog)			
Using pressure	0~1.2kg/cm ²			
Wheel	Silicon carrier w/ fixed screw frame			
Electronic power	AC 220v, 60Hz Single-phase			
Safety device	Safely pressure valve, over temp. controller N.F.B. Buzzer,switch,-etc-			



PT PURNAMA LABORATORY



Pesona Angrek Harapan Blok E.4 No.1
Bekasi Utara 17124

Phone : 021-8838 5761 | Fax : 021-8838 5759

Mobile : 0821 22 555 339 | Email : admin@purnamalab.com

Website : www.purnamalab.com